Cleanser A

Reinigungslösung A
Solución de limpieza A

Order information:

Catalog no.	Content
ULS2100	1x 500 mL
ULS2101	1x 200 mL

Intended use:

Cleanser A is a cleaning solution for use on the Urilyzer® Sed Microscopic Urine Sediment Analyzer. After each measurement, the instrument is cleaned using Cleanser A.

Concentrate ingredients:

Detergent	< 1 %
Buffer	> 5 %
Stabilizers	< 1 %
Isopropanol	> 60%

Stability:

Unopened bottle:	Up to the expiration date	at 2 – 35°C
Opened bottle:	2 months	at 15 – 30°C

Do not freeze! Keep protected from light!

Preparation:

The cleaning solution is concentrated. Dilute 1:10 with distilled water prior to use. Precipitate or floc in the solution is a normal phenomenon.

#ULS2100 (500 mL concentrate):

1. Fill 2.5L distilled water into a 5L container.
2. Slowly pour the concentrated solution into the container.
3. Use 2L distilled water to wash the concentrated solution bottle thoroughly (ensure that no solid residual remains in the bottle), and empty the bottle into the container.
4. Mix and dissolve solid residual.
5. Connect the Cleanser A container to the Urilyzer® Sed with a silicone tube. Be sure to use the correct inlet.

#ULS2101 (200 mL concentrate):

1. Fill 1L distilled water into a 5L container.
2. Slowly pour the concentrated solution into the container.
3. Use 0.8L distilled water to wash the concentrated solution bottle thoroughly (ensure that no solid residual remains in the bottle), and empty the bottle into the container.
4. Mix and dissolve solid residual.
5. Connect the Cleanser A container to the Urilyzer® Sed with a silicone tube. Be sure to use the correct inlet.

Notes:

For in vitro diagnostic use. For professional use only.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal:

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

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