# Cleanser B

Reinigungslösung B  
Solución de limpieza B

## Order information:

Catalog no.	Content
ULS2200	4 x 40 mL
ULS2201	1 x 20 mL

## Intended use:

Cleanser B is a cleaning solution for use on the Urilyzer<sup>®</sup> Sed Microscopic Urine Sediment Analyzer. Cleanser B is used for intensive cleaning and for cleaning the instrument before shutdown.

## Concentrate ingredients:

Detergent	< 5 %
Stabilizers	< 1 %
Sodium Hypochloride	> 60%
Isopropanol	< 5%

## Stability:

Unopened bottle:	Up to the expiration date	at 2 – 35°C
Opened bottle:	2 months	at 15 – 30°C

**Do not freeze! Keep protected from light!**

## Preparation:

The cleaning solution is concentrated. Dilute 1:10 with distilled water prior to use. Precipitate or floc in the reagent is a normal phenomenon.

### #ULS2200 (4 x 40 mL):

1. Fill 200mL distilled water into a 500mL container.
2. Slowly pour the concentrated reagent into the container.
3. Use 160mL distilled water to wash the concentrated reagent bottle thoroughly (ensure that no solid residual remains in the bottle) and empty the bottle into the container.
4. Mix and dissolve solid residual.
5. Connect the Cleanser B container to the Urilyzer<sup>®</sup> Sed with a silicone tube. Be sure to use the correct inlet.

### #ULS2201 (1 x 20 mL):

1. Fill 80mL distilled water into a 500mL container.
2. Slowly pour the concentrated reagent into the container.
3. Use 100mL distilled water to wash the concentrated reagent bottle thoroughly (ensure that no solid residual remains in the bottle) and empty the bottle into the container.
4. Mix and dissolve solid residual.
5. Connect the Cleanser B container to the Urilyzer<sup>®</sup> Sed with a silicone tube. Be sure to use the correct inlet.

## Notes:

For in vitro diagnostic use. For professional use only.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Exercise the normal precautions required for handling all laboratory reagents.

## Disposal:

Please note the legal regulations.

Text passages with grey background were changed in the latest revision of this package insert.

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