UREA HL

UREA HIGH LINEARITY

Order information:

Catalog No.	Contents					
4905	R1	4 x	50ml	R2	4 x	for 50 ml
	R4	1 x	5 ml			
4901	R1	4 x	100ml	R2	4 x	for 100 ml
	R4	1 x	5 ml			

Intended use:

Kinetic in vitro test for the quantitative determination of urea in human serum and plasma.

Summary:



The determination of urea is the most widely used test for the evaluation of kidney function. The test is frequently used in conjunction with the determination of creatinine for the differential diagnosis of prerenal hyperuremia (cardiac decompensation, water depletion increased protein catabolism), renal hyperuremia (glomerulonephritis, chronic nephritis, polycystic kidney, nephrosclerosis, tubular necrosis) and postrenal hyperuremia (obstructions of the urinary tract).

Urea is the final degradation product of protein and amino acid metabolism. In protein catabolism the proteins are broken down to amino acids and deaminated. The ammonia formed in this process is synthesized to urea in the liver. This is the most important catabolic pathway for eliminating excess nitrogen in the human body.

In 1914 Marshall introduced an assay based on the enzyme urease for determining urea in blood. The ammonia released from urea by urease was measured titrimetrically. Numerous other techniques have since been employed to measure the ammonia produced. These include Bertholot's indophenol assay and the reaction of ammonia with Nessler's reagent. Subsequent modifications have been published by Fawcett and Scott and by Chaney and Marbach. In 1965, Talke and Schubert published a totally enzymatic procedure for the determination of urea using the coupled urease/glutamate dehydrogenase (GLDH) enzyme system.

Test principle:

Urea is hydrolysed in presence of urease to form ammonia and CO₂. The ammonia produced reacts with 2 - ketoglutarate and NADH in presence of GLDH to yield glutamate and NAD.

The decrease in absorbance due to the decrease of NADH concentration is proportional to the Urea concentration.

Reagent Concentration:

R1:	
Tris buffer pH 8.0	80 mmol/l
2-Oxoglutarate	15 mmol/l
R2:	
GLDH	6000 U/l
NADH	0.32 mmol/l
Urease	> 100 U/ml
R4:	
Urea	50 mg/dl

Preparation and stability:

Unopened vials are stable until the end of the expiration date when stored at 2–8°C. Dissolve enzyme reagent/R2 with buffer /R1 by gently moving for 20 min. The working solution is stable up to:

5 days	at +20°C to +25°C
28 days	at +2°C to +8°C

Specimen:

Collect serum using standard sampling tubes. Li-heparin, Na-heparin or K-EDTA-Plasma (do not use ammonium heparin). Stability:

7 days	at +20°C to +25°C
7 days	at +2°C to +8°C
1 year	at –20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.

Warning! R2 contains hazardous components.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ± 10% of initial values

Serum/plasma

Icterus: No significant interference up to an index I of 13 (approximate conjugated and unconjugated bilirubin concentration 13 mg/dl)

Hemolysis: No significant interference up to an Index H of 1100 (approximate hemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): No significant interference up to an index L of 225 (approximate triglycerides concentration: 450 mg/dl).

There is poor correlation between turbidity and triglycerides concentration.

Ammonia produced on the cuvette during a GLDH or lactate UV determination interferes with the urea/BUN assay. The urea/BUN must therefore not be installed on the analyzers together with reagents for the GLDH or lactate UV test. In urine endogenous ammonium ions interfere with the urea/BUN assay. Elevated concentrations can occur under acidic conditions (e.g. acidosis.)

Great care must be taken to prevent ammonia contamination of the specimens and calibrators to be analyzed for urea/urea nitrogen.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

• Working solutions as described above

Additional materials required

- Controls as indicated below
- 0.9% NaCl

Manual procedure:

Wavelength:	340 nm, Hg 334 nm
Temperature:	+25 / +37°C
Cuvette:	1 cm light path
Zero adjustment:	air or distilled water

+25°C Semi-Micro

Sample or calib./stand.	40 µl
Working sol.	1000 µl

+37°C Semi-Micro

Sample or calib./stand.	20 µl
Working sol.	1000 µl

Mix and read A₁ after exactly 90 sec. and read A₂ after additional 120 sec.

Calculation:

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ calib./stand.}} \times \text{Calib./stand. conc.} = \text{conc. Urea in mg/dl}$$

$$\text{mg/dl Urea} \times 0,467 = \text{mg/dl Urea-N}$$

Measuring /reportable range:

Up to 400 mg/dl (66.6 mmol/l)

Expected values:

Serum/Plasma

10 - 50 mg/dl (1.7 - 8.3 mmol/l) Urea

6 - 20 mg/dl Urea - Nitrogen

Reference ranges for children are given in the brochure "Reference ranges for adults and children; preanalytical considerations" by Heil W. Koberstein R, Zawta B. (published by Boehringer Mannheim GmbH1997).

24 hour urine

10 - 35 g/24h (170 - 580 mmol/l/24h) corresponding to 670 - 2300 mg/dl*

(110 - 390 mmol/l)* Urea

12 - 20/24h g Urea - N corresponding to 800-1666 mg/dl

* Based on average urine output of 1.5/24h

The expected values are influenced by daily intake of proteins on relation to the body weight.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the urea results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit):

Detection limit: 5 mg/dl or 0.83 mmol/l

The lower detection limit represents the lowest measurable urea concentration that can be distinguished from zero.

UREA HL

UREA HIGH LINEARITY

Imprecision:

serum

Reproducibility was determined using controls (n = 20). The following results were obtained:

within run			
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	42.53	1.78	4.19
Contronorm	103.80	2.88	2.77
Controptath	150.40	3.80	2.53

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	41.98	1.31	3.12
Contronorm	150.68	5.61	3.72
Controptath	554.92	27.12	4.89

Method comparison:

A comparison of the Analyticon UREA-HL (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1.009x - 0.144; \quad r = 0.999$$

Quality Control:

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal® E	10 x 3 ml	#1430
Bio Cal®	20 x 3 ml	#1420

R4: calibrator provided in kit

Disposal:

Please note the legal regulations.

Literature:

- Berthelot MPE Repert Chim Appl 1859;282
- Chaney AL Marbach EP Clin Chem 1962;8:131
- Colombo J-P (ed) Klinisch - chemische Urindiagnostik. Rotkreuz: Labolive-Verlagsgesellschaft 1994:180
- Fawcett JK Scott JE J Clin Chem 1962;13:156
- Gentzkow CJ. J Biol Chem 1942;143:531
- Glick MR. Ryder KW. Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation, Clin Chem 1986;32:470-474
- Guder WG, Narayanan S, Wisser H, Zawta B, List of Analytes ; Pre analytical Variables. Brochure in Samples: From the Patient to the Laboratory. Darmstadt :GIT Verlag 1996
- Krieg M et al J Clin Chem Clin Biochem 1986;24:863
- Mac Kay EM Mac Kay LL J clin Invest 1927;4:295
- Marshall EK Jr. Biol Chem 1913;15:487
- Neumann U, Ziegenhorn J, Scand J Clin Lab Invest 1977;37, Supplement 147: Abstract 97
- Talke H Schubert GE Enzymatische Harnstoffbestimmung im Blut und Serum im optischen Test nach Warburg. Klin Wschr 1965;43:174-175
- Tietz NW. Clinical Guide to Laboratory Tests 3rd ed. Philadelphia, Pa. WB Saunders Company. 1995:624
- Tietz NW. Clinical Guide to Laboratory Tests 3rd ed. Philadelphia, Pa. WB Saunders Company. 1995:622 - 629
- Tietz NW. Fundamentals of Clinical Chemistry Philadelphia, Pa: WB Saunders Co 1976:991

Text passages with grey background were changed in the latest revision of this package insert.

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