Order information:

Catalog No.		Contents			
H9705	Hit I (ILab*)	R1	6 x 40 ml	R2	3 x 49 ml
H9703	Hit 917 (AU*)	R1	6 x 60 ml	R2	6 x 40 ml
AU9703	AU	R1	6 x 60 ml	R2	6 x 40 ml

(*) Kit contains only reagent barcode for Hitachi systems.

System information:

Hitachi 911/917: ACN 418 (Serum/Plasma)
Hitachi 917: ACN 470 (Urine)

For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

Intended use:

Enzymatic *in vitro* test for the quantitative determination of urea in human serum, plasma and urine.



Summary:

The determination of urea is the most widely used test for the evaluation of kidney function. The test is frequently used in conjunction with the determination of creatinine for the differential diagnosis of prerenal hyperuremia (cardiac decompensation, water depletion increased protein catabolism), renal hyperuremia (glomerulonephritis, chronic nephritis, polycystic kidney, nephrosclerosis, tubular necrosis) and postrenal hyperuremia (obstructions of the urinary tract).

Urea is the final degradation product of protein and amino acid metabolism. In protein catabolism the proteins are broken down to amino acids and deaminated. The ammonia formed in this process is synthesized to urea in the liver. This is the most important catabolic pathway for eliminating excess nitrogen in the human body. In 1914 Marshall introduced an assay based on the enzyme urease for determining urea in blood. The ammonia released from urea by urease was measured titrimetrically. Numerous other techniques have since been employed to measure the ammonia produced. These include Bertholot's indophenol assay and the reaction of ammonia with Nessler's reagent. Subsequent modifications have been published by Fawcett and Scott and by Chaney and Marbach. In 1995, Talke and Schubert published a totally enzymatic procedure for the determination of urea using the coupled urease/glutamate dehydrogenase (GLDH) enzyme system. The Analyticon UREA/BUN assay is based on the completely enzymatic method. It has been optimized for automatic analyzers that permit kinetic measurements.

Test principle:

Urea is hydrolysed in presence of urease to produce ammonia and CO₂. The ammonia produced combines with 2-oxoglutarate and NADH in presence of GLDH to yield glutamate and NAD.

The decrease in absorbance due to consumption of NADH is measured kinetically.

Reagent concentration:

R1:
BICIN* buffer pH 7.6 50 mmol/l
GLDH ≥ 0.80 U/l
Urease ≥ 12 U/ml

R2:
TRIS** buffer pH 9.6 100 mmol/l
2-oxoglutarate 8.3 mmol/l
NADH ≥ 0.23 mmol/l

* BICIN = N,N-bis(2-hydroxyethyl)-glycine

** TRIS = Tris (hydroxymethyl)-aminomethane

Preparation and stability:

R1: ready for use

R2: ready for use

The reagents are stable up to the expiry date on the label when stored at +2°C to +8°C.

Onboard stability: R1 28 days
R2 28 days

Notes:

For *in vitro* diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Specimen:

Serum

Collect serum using standard sampling tubes.

Li-heparin, Na-heparin or K-EDTA plasma. Do not use ammonium heparin.

Stability: 7 days at +20°C to +25°C
7 days at +2°C to +8°C
1 year at -20°C

Urine

Collect urine without using preservatives.

Stability: 2 days at +20°C to +25°C
7 days at +2°C to +8°C
1 month at -20°C

- Urine samples are automatically diluted 1 + 19 with 0.9% NaCl or distilled water in the analyzer. The respective dilutions are taken into account in the calculation of results.

- Analyzers without automatic sample dilution

Manually dilute urine samples with 0.9% NaCl or distilled water

(e.g. 1 + 10). Multiply the result by the appropriate dilution factor (e.g. 11).

Centrifuge samples containing precipitate before performing the assay.

Limitations - interference:

Criterion: Recovery within ±10% of initial value.

Icterus: No significant interference up to an I index of 100 (approximate conjugated and unconjugated bilirubin concentration 100 mg/dl)

Hemolysis: No significant interference up to an H Index of 800 (approximate hemoglobin concentration: 800 mg/dl).

Lipemia (Intralipid): No significant interference up to an L index of 1200 (approximate triglycerides concentration: 2400 mg/dl).

There is poor correlation between turbidity and triglycerides concentration

Ammonia produced on the cuvette during a GLDH or lactate UV determination interferes with the UREA/BUN assay. The urea/BUN must therefore not be installed on the analyzers together with reagents for the GLDH or lactate UV test. In urine endogenous ammonium ions interfere with the urea/BUN assay. Elevated concentrations can occur under acidic conditions (e.g. acidosis.)

Great care must be taken to prevent ammonia contamination of the specimens and calibrators to be analyzed for urea/urea nitrogen.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:

For this presentation is no manual procedure available. If the use of a manual procedure is necessary, please use Fluitest UREA UV Art, #4705 or #4701.

Calculation:

Conversion into SI-units, relation between urea and urea-nitrogen:

$$\text{mg/dl} \times 0,166 = \text{mmol/l (urea)}$$

$$\text{mg/dl urea} \times 0,467 = \text{mg/dl urea-nitrogen}$$

Measuring/reportable range:

Serum/plasma

5 - 400 mg/dl (0.83 to 66.4 mmol/l) urea or 2 - 186 mg/dl urea nitrogen.

Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl or distilled water (e.g. 1+2). Multiply the result by the appropriate dilution factor (e.g. 3).

Calculated extended ranges with rerun for 900-series analyzer:

- Roche/Hitachi 911:

Up to 800 mg/dl (133 mmol/l) urea or 374 mg/dl urea nitrogen.

- Roche/Hitachi 917:

Up to 600 mg/dl (99.6 mmol/l) urea or 280 mg/dl urea nitrogen.

Urine

- Roche/Hitachi 911/917:

50 - 8000 mg/dl (8.3-1328 mmol/l) urea or 23-3736mg/dl urea nitrogen.

Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl or distilled water (e.g. 1 + 2). Multiply the result by the appropriate dilution factor (e.g. 3).

Expected values:

Serum/plasma

10 – 50 mg/dl (1.7 – 8.3 mmol/l) urea

Reference ranges for children are given in the brochure "Reference ranges for adults and children; preanalytical considerations" by Heil W. Koberstein R, Zawta B. (published by Boehringer Mannheim GmbH 1997).

Morning urine

847 – 2967 mg/dl (141 – 494 mmol/l)

24-hour urine

10 – 35 g/24h (170 – 580 mmol/24h), corresponding to

670 – 2300 mg/dl (110 – 390 mmol/l)

The expected values are influenced by daily take-up of proteins on relation to the body weight.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the urea results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 5 mg/dl (0.83 mmol/l)

The lower detection limit represents the lowest measurable urea activity that can be distinguished from zero.

Precision:

Reproducibility within run was determined using human samples and controls (n = 20). The following results were obtained:

Serum			
Sample	Within run		
	Mean mg/dl	SD mg/dl	CV %
Sample1	41.8	1.34	3.21
Sample2	98.4	1.63	1.66
Sample3	142	3.03	2.13

Serum			
Sample	Between day		
	Mean mg/dl	SD mg/dl	CV %
Sample1	43.5	1.34	3.07
Sample2	65.6	2.05	3.12
Sample3	140.2	3.26	2.32

Method comparison:

A comparison of the Analyticon Fluitest® UREA/BUN (y) with a commercial obtainable assay (x) gave the following result:

$$y = 0,993 x + 0,389; \quad r = 0,998$$

Quality Control:

Human Control Serum:

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

Human Urine Control:

Urine control Set	8 x 5 ml	#1507
-------------------	----------	-------

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S1: 0.9% NaCl

S2: Bio Cal®	20 x 3 ml	#1420
Bio Cal® E	10 x 3 ml	#1430

Disposal:

Please note the legal regulations.

Literature:

- Berthelot MPE Reper Chim Appl 1859;282
- Chaney AL Marbach EP Clin Chem 1962;8:131
- Colombo J-P (ed) Klinisch - chemische Urindiagnostik. Rotkreuz: Labolive-Verlagsgesellschaft 1994:180
- Fawcett JK Scott JE J Clin Chem 1962;13:156
- Gentzkow CJ. J Biol Chem 1942;143:531
- Glick MR. Ryder KW. Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation, Clin Chem 1986;32:470-474
- Guder WG, Narayanan S, Wisser H, Zawta B, List of Analytes ; Pre analytical Variables. Brochure in Samples: From the Patient to the Laboratory. Darmstadt :GIT Verlag 1996
- Krieg M et al J Clin Chem Clin Biochem 1986;24:863
- Mac Kay EM Mac Kay LL J clin Invest 1927;4:295
- Marshall EK Jr. Biol Chem 1913;15:487
- Neumann U, Ziegenhorn J, Scand J Clin Lab Invest 1977;37, Supplement 147: Abstract 97
- Talke H Schubert GE Enzymatische Harnstoffbestimmung im Blut und Serum im optischen Test nach Warburg. Klin Wschr 1965;43:174-175
- Tietz NW. Clinical Guide to Laboratory Tests 3rd ed. Philadelphia, Pa. WB Saunders Company. 1995:624, 622 - 629
- Tietz NW. Fundamentals of Clinical Chemistry Philadelphia, Pa: WB Saunders Co 1976:991

Text passages with grey background were changed in the latest revision of this package insert.

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
* BICIN = N,N-Bis(2-hydroxyethyl)-Glycin

** TRIS = Tris (hydroxymethyl)-Aminomethan

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