Order information

| Catalog No. | Contents | | | | | |
|-------------|----------|-----|--------|----|-----|-----------------|
| 5113 | R1 | 1 x | 100 ml | R2 | 1 x | 1 ml for 100 ml |
| | R3 | 1 x | 25 ml | R4 | 1 x | 5 ml |
| 5111 | R1 | 3 x | 100 ml | R2 | 3 x | 1 ml for 100ml |
| | R3 | 1 x | 80 ml | R4 | 1 x | 5 ml |

Intended use:

Enzymatic in vitro test for the quantitative determination of urea in human serum, plasma and urine.

Summary:

The determination of urea is the most widely used test for the evaluation of kidney function. The test is frequently used in conjunction with the determination of creatinine for the differential diagnosis of prerenal hyperuremia (cardiac decompensation, water depletion increased protein catabolism), renal hyperuremia (glomerulonephritis, chronic nephritis, polycystic kidney, nephrosclerosis, tubular necrosis) and postrenal hyperuremia (obstructions of the urinary tract).

Urea is the final degradation product of protein and amino acid metabolism. In protein catabolism the proteins are broken down to amino acids and deaminated. The ammonia formed in this process is metabolized to urea in the liver. This is the most important catabolic pathway for eliminating excess nitrogen in the human body. In 1914 Marshall introduced an assay based on the enzyme urease for determining urea in blood. The ammonia released from urea by urease was measured titrimetrically. Numerous other techniques have since been employed to measure the ammonia produced. These include Berthelot's indophenol assay and the reaction of ammonia with Nessler's reagent. Subsequent modifications have been published by Fawcett and Scott and by Chaney and Marbach. In 1995, Talke and Schubert published a totally enzymatic procedure for the determination of urea using the coupled urease/glutamate dehydrogenase (GLDH) enzyme system. The Analyticon UREA, col. assay is based on Berthelot method.

Test principle:

Enzymatic determination of Urea according to the following reaction:

The ammonium ions formed react with salicylate and hypochloride to give a green dye (2,2-Dicarboxylindophenol).

Reagent concentration:

| | |
|-------------------------|-------------------------|
| R1: | |
| Phosphate buffer pH 6.7 | 60 mmol/l |
| EDTA | 1.5 mmol/l |
| Sodium salicylate | 60 mmol/l |
| Sodium nitroprusside | 5.2 mmol/l |
| R2: | |
| Urease | > 5000 U/l |
| R3: | |
| Sodium hypochloride | 18 mmol/l |
| Sodium hydroxide | 450 mmol/l |
| R4: | |
| Urea | 50 mg/dl (8.325 mmol/l) |

Preparation and stability:

Working reagent :

Add 1 vial R2 to 1 bottle of R1. The working solution is stable up to:

| | |
|---------|------------------|
| 4 weeks | at +2 to +8°C |
| 6 days | at +20 to +25°C. |

R3: ready for use

R4: ready for use

The reagents are stable up to the expiry date specified when stored at +2°C to +8°C.

Protect all reagents against direct light.

For photometers, which need a minimum of 2 ml, dilute R3 by 4 volumes of bidistilled water, e.g. 10 ml of R3 and 40 ml of bidistilled water. The diluted reagent is stable for:

6 months at +2 to +8°C.

Specimen:

Centrifuge samples containing precipitate before performing the assay.

Serum, plasma

Collect serum using standard sampling tubes.

Heparin, citrated plasma

| | | |
|------------|--------|-------------------|
| Stability: | 7 days | at +20°C to +25°C |
| | 7 days | at +2°C to +8°C |
| | 1 year | at -20°C |

Urine

Urine diluted 1:10 with distilled water.

Collect urine without using preservatives.

| | | |
|------------|---------|-------------------|
| Stability: | 2 days | at +20°C to +25°C |
| | 7 days | at +2°C to +8°C |
| | 1 month | at -20°C |

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Criterion: Recovery within ±10% of initial values

Icterus: No significant interference up to an I index of 30 (approximate conjugated and unconjugated bilirubin concentration 30 mg/dl)

Hemolysis: No significant interference up to an H Index of 1100 (approximate hemoglobin concentration: 1100 mg/dl).

Lipemia (Intralipid): Triglycerides may interfere. There is poor correlation between turbidity and triglycerides concentration.

Ammonia produced on the cuvette during a GLDH or lactate UV determination interferes with the urea assay. The urea must therefore not be installed on the analyzers together with reagents for the GLDH or lactate UV test. In urine endogenous ammonium ions interfere with the urea assay. Elevated concentrations can occur under acidic conditions (e.g. acidosis.)

Great care must be taken to prevent ammonia contamination of the specimens and calibrators to be analyzed for urea/urea nitrogen.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below

- 0.9% NaCl

Manual procedure:

| | |
|------------------|-----------------------|
| Wavelength: | 580- 600 nm |
| Temperature: | +20 / +25 / +37°C |
| Cuvette: | 1 cm light path |
| Zero adjustment: | against reagent blank |

| | Reagent blank | Standard | Sample |
|--------------|---------------|----------|---------|
| Working sol. | 1000 µl | 1000 µl | 1000 µl |
| R4/Standard | --- | 10 µl | --- |
| Sample | --- | --- | 10 µl |

Mix, incubate at +37°C for 5 min. or for 10 min at +20 / +25°C.
Then add:

| | | | |
|-----------------------------|---------|---------|---------|
| R3 | 200µl | 200µl | 200µl |
| or diluted reagent/R3 (1:5) | 1000 µl | 1000 µl | 1000 µl |

Mix, incubate at +37°C for 5 min. or for 10 min at +20 / +25 °C and read absorbance against reagent blank (ΔA).

Calculation:

$$\frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Standard}}} \times \text{Standard conc.} = \text{Urea conc. (mg/dl)}$$

Measuring /reportable range:

Up to 200 mg/dl (33.3 mmol/l) for the normal pipetting scheme. The linearity using the diluted reagent/R3 is enhanced up to 300 mg/dl (50 mmol/l)

In case of higher results, dilute sample 1 + 1 with saline solution and repeat test. Multiply result by 2.

Expected values:

Serum

15 – 45mg/dl (1.7 – 8.3 mmol/l)

4.7 – 23 mg/100 ml Urea - N

Reference ranges for children are given in the brochure "Reference ranges for adults and children; preanalytical considerations" by Heil W. Koberstein R, Zawta B. (published by Boehringer Mannheim GmbH 1997).

24 hour urine

20 – 36 g/24h (333 – 600 mmol/l/24h)

9.4 – 16.9 g Urea – N

The expected values are influenced by daily take – up of proteins on relation to the body weight.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the urea results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 1 mg/dl

The lower detection limit represents the lowest measurable urea activity that can be distinguished from zero.

Imprecision:

Serum

Reproducibility was determined using human samples and controls between day (n = 20). The following results were obtained:

| Within run | | | |
|-------------------|---------------|-------------|---------|
| Sample | Mean mg/dl | SD mg/dl | CV % |
| Sample 1 | 39.61 | 0.88 | 2.21 |
| Sample 2 | 90.95 | 5.65 | 6.21 |
| Sample 3 | 139.78 | 1.95 | 1.40 |

Reproducibility was determined using human samples and controls between day (n = 18). The following results were obtained:

| Between day | | | |
|--------------------|---------------|-------------|---------|
| Sample | Mean mg/dl | SD mg/dl | CV % |
| Sample 1 | 39.85 | 1.52 | 3.81 |
| Sample 2 | 89.48 | 3.47 | 3.87 |
| Sample 3 | 140.20 | 5.27 | 3.76 |

Method comparison:

A comparison of the Analyticon Fluitest UREA col. (y) with a commercial obtainable assay (x) gave the following result:

$$y = 1.000 x + 0.423 ; \quad r = 0.998$$

Quality Control:

Human Control Serum:

| | | |
|-------------------|-----------|-------|
| Contronorm® Plus | 5 x 5 ml | #1205 |
| | 20 x 5 ml | #1220 |
| Controptath® Plus | 5 x 5 ml | #1305 |
| | 20 x 5 ml | #1320 |

Human Urine Control:

| | | |
|-------------------|----------|-------|
| Urine control Set | 8 x 5 ml | #1507 |
|-------------------|----------|-------|

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

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