Order information:

Catalog No.	Contents					
4705	R1	4 x	50 ml	R2	4 x	10 ml
	R4	1 x	5 ml			

Intended use:

Enzymatic in vitro test for the quantitative determination of urea in human serum, plasma and urine.


Summary:

The determination of urea is the most widely used test for the evaluation of kidney function. The test is frequently used in conjunction with the determination of creatinine for the differential diagnosis of prerenal hyperuremia (cardiac decompensation, water depletion increased protein catabolism), renal hyperuremia (glomerulonephritis, chronic nephritis, polycystic kidney, nephrosclerosis, tubular necrosis) and postrenal hyperuremia (obstructions of the urinary tract). Urea is the final degradation product of protein and amino acid metabolism. In protein catabolism the proteins are broken down to amino acids and deaminated. The ammonia formed in this process is synthesized to urea in the liver. This is the most important catabolic pathway for eliminating excess nitrogen in the human body.

In 1914 Marshall introduced an assay based on the enzyme urease for determining urea in blood. The ammonia released from urea by urease was measured titrimetrically. Numerous other techniques have since been employed to measure the ammonia produced. These include Bertholot's indophenol assay and the reaction of ammonia with Nessler's reagent. Subsequent modifications have been published by Fawcett and Scott and by Chaney and Marbach. In 1995, Talke and Schubert published a totally enzymatic procedure for the determination of urea using the coupled urease/glutamate dehydrogenase (GLDH) enzyme system. The Analyticon UREA-UV assay is based on the completely enzymatic method. It has been optimized for automatic analyzers that permit kinetic measurements.

Test principle:

Urea is hydrolysed in presence of urease to produce ammonia and CO₂. The ammonia produced combines with 2-oxoglutarate and NADH in presence of GLDH to yield glutamate and NAD.

The decrease in absorbance due to the decrease of NADH concentration in unit time is proportional to the urea concentration.

Reagent concentration:

R1:	
Tris buffer pH 8.0	50 mmol/l
2-oxoglutarate	15 mmol/l
Urease	≥ 1000 U/l
GLDH	≥ 6000 U/l
R2:	
NADH	0.18 mmol/l
R4:	
Urea	50 mg/dl (8.325 mmol/l)

Preparation and stability:

R1: ready for use
 R2: ready for use
 R4: ready for use
 The reagents are stable up to the expiry date on the label when stored at 2-8°C.
 Onboard stability: R1 28 days
 R2 28 days

Alternatively the reagent can be used as single reagent (manual procedure). Mix 5 volumes enzyme reagent/R1 with one volume starting buffer /R2. The working solution is stable up to:

14 day	at + 2°C to +8°C
3 days	at +20°C to +25°C

Specimen:

Serum
 Collect serum using standard sampling tubes. Heparin, citrated plasma.
 Stability: 7 days at +20°C to +25°C
 7 days at +2°C to +8°C
 1 year at -20°C

Urine
 Urine diluted 1:50 with distilled water.
 Collect urine without using preservatives.
 Stability: 2 days at +20°C to +25°C
 7 days at +2°C to +8°C
 1 month at -20°C

Centrifuge samples containing precipitate before performing the assay.

Notes:

For in vitro diagnostic use.
 The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Limitations - interference:

Ammonia produced on the cuvette during a GLDH or lactate UV determination interferes with the urea-UV assay. The urea-UV must therefore not be installed on the analyzers together with reagents for the GLDH or lactate UV test. In urine endogenous ammonium ions interfere with the urea-UV assay. Elevated concentrations can occur under acidic conditions (e.g. acidosis.)

Great care must be taken to prevent ammonia contamination of the specimens and calibrators to be analyzed for urea/urea nitrogen.

Blood samples should only be drawn prior to the administration of sulfasalazine and sulfapyridine. Physiological plasma concentrations of sulfasalazine or sulfapyridine can lead to incorrect results.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure:

- Materials provided**
- Working solutions as described above
- Additional materials required**
- Calibrators and controls as indicated below
 - 0.9% NaCl

Manual procedure		
Wavelength:	340nm, Hg 334 nm	
Temperature:	+25°C or +37°C	
Cuvette:	1 cm light path	
Zero adjustment:	distilled water	
	Standard/Calibrator	Sample
Working reagent	1000 µl	1000 µl
Standard R4	10 µl	----
Sample	----	10 µl
Mix, and read A ₁ after exactly 30 sec., read A ₂ after additional 60 sec.		
Calculation:		
Determine the absorbance change as		
ΔA sample = [A ₁ (sample) – A ₂ (sample)]		
ΔA calibrator = [A ₁ (calibrator)– A ₂ (calibrator)]		
and use this for calculation.		
$\frac{\Delta A \text{ sample}}{\Delta A \text{ calibrator}} \times \text{calibrator conc.} = \text{Urea conc. (mg/dl)}$		

Measuring/reportable range

Kinetic measuring: linear up to 400 mg/dl (66.7 mmol/l)
 Conversion in SI-units, relation between urea and urea-nitrogen:
 mg/dl x 0.166 = mmol/l (urea)
 mg/dl urea x 0.467 = mg/dl urea-nitrogen

Expected values:

Serum
 15 – 45 mg/dl (1.7 – 8.3 mmol/l) urea
 4.7 – 2300 mg/dl (31.1-15.3 mmol/l) urea-nitrogen
 Reference ranges for children are given in the brochure "Reference ranges for adults and children; preanalytical considerations" by Heil W. Koberstein R, Zawta B. (published by Boehringer Mannheim GmbH 1997).

24-hour urine
 20 – 36 g/l (333 – 600 mmol/l) urea
 9.4 – 16.9 g/l (6.3-11.3 mmol/l) urea-nitrogen

The expected values are influenced by daily take-up of proteins on relation to the body weight.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes the urea results should always be assayed in conjunction with the patient's medical history, clinical examinations and other findings.

Analytical sensitivity (lower detection limit)

Detection limit: 0.5 mg/dl
 The lower detection limit represents the lowest measurable urea activity that can be distinguished from zero.

Imprecision:

Urea

Reproducibility was determined using human samples and controls between day. The following results were obtained:

Serum

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	31	1.1	3.4
Contronorm	50	0.9	1.8
Controptath	145	1.6	1.1

Urine

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	692	13.1	1.9
Contronorm	1753	25.6	1.5
Controptath	1813	23.9	1.3

Urea-N

Reproducibility was determined using human samples and controls between day. The following results were obtained:

Serum

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	17	0.56	3.3
Contronorm	20	0.50	2.4
Controptath	66	0.82	1.2

Urine

Between day			
Sample	Mean mg/dl	SD mg/dl	CV %
Human serum	568	17.94	3.2
Contronorm	500	3.68	0.7
Controptath	769	10.58	1.4

Quality Control:

Human Control Serum:

Contronorm® Plus 5 x 5 ml #1205
 20 x 5 ml #1220

Controptath® Plus 5 x 5 ml #1305
 20 x 5 ml #1320

Human Urine Control:

Urine control Set 8 x 5 ml #1507

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

S2: Bio Cal® E 10 x 3 ml #1430
 Bio Cal® 20 x 3 ml #1420

R4: calibrator provided in kit

Disposal:

Please note the legal regulations.

Literature:

- Berthelot MPE Repert Chim Appl 1859;282
- Chaney AL Marbach EP Clin Chem 1962;8:131
- Colombo J-P (ed) Klinisch - chemische Urindiagnostik. Rotkreuz: Labolive-Verlagsgesellschaft 1994;180
- Fawcett JK Scott JE J Clin Chem 1962;13:156
- Gentzkow CJ. J Biol Chem 1942;143:531
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation, Clin Chem 1986;32:470-474
- Guder WG, Narayanan S, Wisser H, Zawta B, List of Analytes ; Pre analytical Variables. Brochure in Samples: From the Patient to the Laboratory. Darmstadt :GIT Verlag 1996
- Krieg M et al J Clin Chem Clin Biochem 1986;24:863
- Mac Kay EM Mac Kay LL J clin Invest 1927;4:295
- Marshall EK Jr. Biol Chem 1913;15:487
- Neumann U, Ziegenhorn J, Scand J Clin Lab Invest 1977;37, Supplement 147: Abstract 97
- Talke H Schubert GE Enzymatische Harnstoffbestimmung im Blut und Serum im optischen Test nach Warburg. Klin Wschr 1965;43:174-175
- Tietz NW. Clinical Guide to Laboratory Tests 3rd ed. Philadelphia, Pa. WB Saunders Company. 1995:624, 622 - 629
- Tietz NW. Fundamentals of Clinical Chemistry Philadelphia, Pa: WB Saunders Co 1976:991

Text passages with grey background were changed in the latest revision of this package insert.

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