BioLyzer® Order Information:

Catalog No.	BioLyzer	Contents
B9301	200 / 600	R1 6 x 50 ml
B9303	300 / 600*	R1 6 x 60 ml

*only for instruments with a reagent tray for 20 / 70ml bottles

Intended use:

For the quantitative determination of protein in urine and cerebrospinal fluid on automated clinical chemistry analyzers.

Summary:

Urine is formed by ultrafiltration of plasma across the glomerular capillary wall. Proteins with a relative molecule mass > 40 000 are almost completely retained, while smaller substances easily enter the glomerular filtrate.

Most CSF protein originates by diffusion from plasma across the blood-CSF barrier. Elevated levels occur as a result of increased permeability of the blood-CFS barrier or with increased local synthesis of immunoglobulins.

Turbidimetric methods using trichloroacetic acid (TCA) or sulfosalicylic acid (SSA) require precipitation of the protein in the sample; the resulting turbidity may be unstable and flocculate.

Protein measurements in urine are used in the diagnosis and treatment of disease conditions such as renal or heart diseases, or thyroid disorders, which are characterized by proteinuria or albuminuria.

CSF protein measurements are used in the diagnosis and treatment of conditions such as meningitis, brain tumors, and infections of the central nervous systems.

Test principle:

Pyrogallol-red-molybdate complexes are bound to proteins. The binding causes an increase of absorbance at 600 nm. This increase is direct proportional to the protein concentration. Due to the presence of detergents in the reagents, different types of proteins give similar recoveries.

Reagent Concentration:

R1:	
Pyrogallol red	0.035 mmol/l
Sodium molybdate	0.018 mmol/l
Succinate buffer	50 mmol/l
Detergent	2 %

Preparation and stability:

Reagent and standard are ready to use.

The reagents are stable up to the stated expiry date when stored at +2 to +8°C and protected from light.

On board stability: 28 days
Contamination after opening must be avoided! Protect from light.

Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Specimen:

Centrifuge samples containing precipitate before performing the assay.

Collection and preparation

Urine:

Use random or 24 hour urine specimens. Use no preservatives.

Refrigerate specimen during collection.

Cerebrospinal Fluid (CSF):

No special additives are required. Blood in a CSF specimen invalidates the protein value.

Specimens may be stored at temperatures between +2°C and +8°C for 48 hours.

Limitations – interference:

Criterion: recovery ± 10 % of initial values

Icterus: No significant interference up to an index I of 15 (approximate conjugated and unconjugated bilirubin concentration: (15 mg/dl).

Hemolysis: Due to protein nature of haemoglobin hemolysis leads elevated values dependent on the degree of lysis of the erythrocytes.

No significant interference from:

Ascorbic Acid	Creatinine	Glucose
Phosphorus	Urea	Magnesium
Sodium Citrate	Caffeine	Cefazolin Sodium
Chlorpromazine	Calcium	L-Dopa Gentamicin
Sulfate	Sodium Oxalate	Uric Acid

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

Testing procedure

Materials provided

• Working solutions as described above

Additional materials required

• Calibrators and controls as indicated below
• 0.9% NaCl

Measuring /reportable range:

5 – 100 mg/dl (0.05 – 1 g/l)

Determine samples having higher concentrations via the rerun function using 0.9 % NaCl.

Analytical sensitivity (lower detection limit):

Detection limit: 5 mg/dl or 0.05 g/l

The lower detection limit represents the lowest measurable USP Protein concentration that can be distinguished from zero. It is calculated as the concentration lying three standard deviations above that of the lowest standard.

Expected values:

Urine:

Random	< 12 mg/dl (< 120 mg/l)
24 hr	< 150 mg/24h (<150 mg/d)

Cerebrospinal Fluid:

Adults: < 45 mg/dl (< 0.45 g/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the USP results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol. The following results were obtained.

Urine	Within run		
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Control 1	45.13	0.96	2.1
Control 2	82.11	1.42	1.7
Control 3	103.62	1.24	1.2

Urine	Run to run		
Sample	Mean (mg/dl)	SD (mg/dl)	CV %
Control 1	72.01	2.80	3.9
Control 2	74.20	2.14	2.9
Control 3	103.21	3.38	3.3

Method comparison:

A comparison of the Analyticon Fluitest USP (y) with a commercial obtainable assay (x) gave with the following result (mg/dl):

$$y = 0.8513x + 1.0211; r = 0.976$$

Quality Control:

Human Urine Control:

Urine Control Set 8 x 5 ml #1507

Calibration:

Calibration Type: Linear

S1: 0.9% NaCl

S2: Protein- C Standard 3 x 5 ml #1491

Calibration frequency:

Calibration is recommended

• After bottle change
• As required following quality control procedures

Disposal:

Please note the legal regulations.

Literature:

- Iwata I, Nishikaze O. Clin Chem. 1979;25/7:1317 – 1319.
- Koumantakis G. Fluorescein Interference with Urinary Creatinine and Protein Measurements. Clin Chem. 1991;37/10:1799.
- Luxton R, Patel P, Keir G, Thompson E. Clin Chem. 1989; 35/8: 1731 – 1734.
- Tietz NW. clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, Pa. WB Saunders co; 1995;520.
- Tietz NW. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia, Pa: WB Saunders co; 1987:336,339,340,341.
- Watanabe N, Kamei S et al., Urinary protein as measured with a pyrogallol red-molybdate complex manually and in a Hitachi 726 automated analyzer. Clin Chem 1986; 32:1551-1554.

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Literatur:

1. Iwata I, Nishikaze O. Clin Chem. 1979;25/7:1317 – 1319.
2. Koumantakis G. Fluorescein Interference with Urinary Creatinine and Protein Measurements. Clin Chem. 1991;37/10:1799.
3. Luxton R, Patel P, Keir G, Thompson E. Clin Chem. 1989; 35/8:1731– 1734.
4. Tietz NW. clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, Pa. WB Saunders co; 1995;520.
5. Tietz NW. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia, Pa: WB Saunders co; 1987:336,339,340,341.
6. Watanabe N, Kamei S et al., Urinary protein as measured with a pyrogallol red-molybdate complex manually and in a Hitachi 726 automated analyzer. Clin Chem 1986; 32:1551-1554.