Order Information:

Catalog No.	Contents		
915	R1 6 x 50 ml	R4 1 x 5 ml	
AU9303 AU	R1 6 x 60 ml		

Intended use:

For the quantitative determination of protein in urine and cerebrospinal fluid on automated clinical chemistry analyzers.

Summary:

Urine is formed by ultrafiltration of plasma across the glomerular capillary wall. Proteins with a relative molecule mass > 40 000 are almost completely retained, while smaller substances easily enter the glomerular filtrate. Most CSF protein originates by diffusion from plasma across the blood-CSF barrier. Elevated levels occur as a result of increased permeability of the blood-CFS barrier or with increased local synthesis of immunoglobulins. Turbidimetric methods using trichloroacetic acid (TCA) or sulfosalicylic acid (SSA) require precipitation of the protein in the sample; the resulting turbidity may be unstable and flocculate. Dye-binding methods such as Protein measurements in urine are used in the diagnosis and treatment of disease conditions such as renal or heart diseases, or thyroid disorders, which are characterized by proteinuria or albuminuria. CSF protein measurements are used in the diagnosis and treatment of conditions such as meningitis, brain tumours, and infections of the central nervous systems.

Test principle:

Pyrogallol-red-molybdate complexes are bound to proteins. The binding causes an increase of absorbance at 600nm. This increase is direct proportional to the protein concentration. Due to the presence of detergents in the reagents, different types of proteins give similar recoveries.

Reagent Concentration:

R1:	
Pyrogallol red	0.035 mmol/l
Sodium molybdate	0.018 mmol/l
Succinate buffer	50 mmol/l
Detergent	2 %
R4 (#915):	
Protein	100 mg/dl (1 g/l)

According to a composition of 70% albumin and 30% immunoglobulins

Preparation and stability:

Reagent and standard are ready to use. The reagents are stable up to the stated expiry date when stored at +2 to +8°C and protected from light. Contamination after opening must be avoided! Protect from light.

Notes:

For in vitro diagnostic use. The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

Specimen:

Centrifuge samples containing precipitate before performing the assay.

Collection and preparation

Urine:
Use random or 24 hour urine specimens. Use no preservatives. Refrigerate specimen during collection.

Cerebrospinal Fluid (CSF):
No special additives are required. Blood in a CSF specimen invalidates the protein value.

Specimens may be stored at temperatures between +2°C and +8°C for 48 hours.

Limitations – interference:

Criterion: recovery ± 10 % of initial values
Icterus: No significant interference up to an index I of 20 (approximate conjugated and unconjugated bilirubin concentration: (20 mg/dl).
Hemolysis: Due to protein nature of haemoglobin hemolysis leads elevated values dependent on the degree of lysis of the erythrocytes.
The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

No significant interference from:

Ascorbic Acid	Creatinine	Glucose
Phosphorus	Urea	Magnesium
Sodium Citrate	Caffeine	Cefazolin Sodium
Chlorpromazine	Calcium	L-Dopa
Sulfate	Sodium Oxalate	Uric Acid
Gentamicin		

Testing procedure

Applications for automated systems are available on request.

Materials provided

- Working solutions as described above
- Additional materials required
- Calibrators and controls as indicated below
- 0.9% NaCl

Manual procedure:

Wavelength:	Hg 600 nm
Temperature:	+25 / +30 / +37°C.
Cuvette:	1 cm light path
Zero adjustment:	Reagent blank; one reagent blank per series only

Macro:

	Blank	Standard	Sample
Standard/ R4	---	75 µl	---
Sample	---	---	75 µl
R1	3000 µl	3000 µl	3000 µl

Micro:

	Blank	Standard	Sample
Standard/ R4	---	25 µl	---
Sample	---	---	25 µl
R1	1000 µl	1000 µl	1000 µl

Mix, incubate for 10 minutes. Measure the absorbance.

Calculation:

$$\frac{\Delta A \text{ sample}}{\Delta A \text{ calibrator}} \times \text{calibrator conc.} = \text{protein conc. in mg/dl}$$

Measuring /reportable range:

5.0 - 160 mg/dl (0.05 - 1.6 g/l) for the micro and macro pipetting schedule. In case of higher concentrations, dilute sample 1 : 4 with saline solution (0.9%) and repeat test. Multiply result by 5.

Analytical sensitivity (lower detection limit)

Detection limit: 2.0 mg/dl or 0.02 g/l
The lower detection limit represents the lowest measurable USP Protein concentration that can be distinguished from zero. It is calculated as the concentration lying three standard deviations above that of the lowest standard.

Expected values:

Urine:	Random	<12 mg/dl (<120 mg/l)
	24 hr	<150 mg/24h (<150 mg/d)
Cerebrospinal Fluid:	Adults:	<45 mg/dl (<0.45 g/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the USP results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Imprecision:

Reproducibility was determined using human samples and controls in an internal protocol. The following results were obtained.

Urine	Between day		
Sample	Mean mg/dl	SD mg/dl	CV %
Human urine	15.11	0.55	3.64
Control 1	19.08	0.69	3.61
Control 2	57.29	0.59	1.03

Urine	Within run		
Sample	Mean mg/dl	SD mg/dl	CV %
Control 1	15.11	0.55	3.64
Control 2	19.08	0.69	3.61

Quality Control:

Human Urine Control: 8 x 5 ml #1507
Urine control Set

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Fluitest[®] USP

ULTRASENSITIVE PROTEIN

Calibration:

S1: 0.9% NaCl

S2: Protein-C Standard 3 x 5 ml #1491
Calibrator R4

Calibration frequency

Recalibration is recommended

- After bottle change
- As required following quality control procedures

Disposal:

Please note the legal regulations.

Literature:

1. Iwata I, Nishikaze O. Clin Chem. 1979;25/7:1317 – 1319.
2. Koumantakis G. Fluorescein Interference with Urinary Creatinine and Protein Measurements. Clin Chem. 1991;37/10:1799.
3. Luxton R, Patel P, Keir G, Thompson E. Clin Chem. 1989; 35/8: 1731 – 1734.
4. Tietz NW. clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, Pa. WB Saunders co; 1995;520.
5. Tietz NW. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia, Pa: WB Saunders co; 1987:336,339,340,341.
6. Watanabe N, Kamei S et al., Urinary protein as measured with a pyrogallol red-molybdate complex manually and in a Hitachi 726 automated analyzer. Clin Chem 1986; 32:1551-1554.

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