# Fluitest® $\alpha$ -HBDH

$\alpha$ -HYDROXYBUTYRATE DEHYDROGENASE

## Order information:

Catalog No.		Contents			
6401	MPR	R1	5 x 20 ml	R2	1 x 25 ml
H6401	Hit I (ILab*)	R1	6 x 20 ml	R2	2 x 13 ml

(\* ) Kit contains only reagent barcode for Hitachi system.

## System information:

Hitachi 911/917: ACN 175  
For application on the Hitachi 912, please use an open system or contact the Analyticon Customer Support.

## Intended use:

In vitro assay for the quantitative determination of  $\alpha$ -hydroxybutyrate dehydrogenase in human serum and plasma.


## Summary:

The  $\alpha$ -hydroxybutyrate dehydrogenase ( $\alpha$ -HBDH) or LDH-1 is an isoenzyme of lactate dehydrogenase (LDH) which uses hydroxybutyrate as additional substrate. Compared to other LDH isoenzymes, the  $\alpha$ -HBDH is predominantly found in the myocardium. Therefore  $\alpha$ -HBDH is used in the diagnostic of myocardial infarction. The quotient HBDH/LDH is used in order to differentiate between kidney and myocardial diseases. A decreased quotient indicates a parenchymal kidney disease while an increased quotient is an indication of the severity and progress of a recent myocardial infarction.

This method is an optimised UV-Test recommended by the German Society for Clinical Chemistry (DGKC) in 1972.

## Test principle:

UV-Test according to DGKC

$\alpha$ -HBDH catalyses the conversion of oxobutyrate to hydroxybutyrate. NADH is oxidized to NAD in the process. The rate of decrease in NADH is directly proportional to the  $\alpha$ -HBDH activity and is measured photometrically.

## Reagent Concentration:

<b>R1:</b>	
Phosphate buffer, pH 7.4	60 mmol/l
2-Oxobutyrate	3.8 mmol/l
<b>R2:</b>	
NADH	>1 mmol/l

## Preparation and stability:

Unopened kit components are stable up to the expiry date when stored at +2°C to +8°C. Protect from light.

### Substrate start:

R1: Ready for use

R2: Ready for use

Onboard stability:	R1	28 days
	R2	28 days

### Sample start:

Gently mix 4 parts R1 with 1 part R2 as required. Protect monoreagent from light.

The working reagent is stable:	8 hours	at +15 to +25°C
	5 days	at +2 to +8°C

## Specimen:

Collect serum using standard sampling tubes.

Serum, heparin- or EDTA-plasma.

Stability:	3 days	at +15 to +25°C
	20 days	at +2 to +8°C

Centrifuge samples containing precipitate before performing the assay.

## Notes:

For in vitro diagnostic use.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>. Exercise the normal precautions required for handling all laboratory reagents.

## Limitations - interference:

Criterion: Recovery within  $\pm 10\%$  of initial value.

Ascorbic acid: No significant interference up to 30 mg/dl.

Icterus: No significant interference up to an index I of 40 (approximate 40 mg/dl bilirubin)

Lipemia (Intralipid): No significant interference up to an index L of 1000 (approximate 2000 mg/dl triglycerides). There is poor correlation between turbidity and triglycerides concentration.

Hemolysis: even small quantities interfere.

The effect of drugs or their metabolic products on the test is not known in all cases. In case of doubt it is recommended to repeat the test after the medication has been discontinued. However, the medication may only be discontinued according to instructions of the attending physician.

## Testing procedure:

Applications for automated systems are available on request.

### Materials provided

• Working solutions as described above

### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

### Manual procedure for sample start:

Wavelength	340 nm (Hg 334 nm)
Temperature	+37°C
Cuvette	1 cm light path
Zero adjustment	Air or distilled water

Serum/Plasma

Working reagent	1000 $\mu$ l
Sample	10 $\mu$ l

Mix, incubate for 60 sec., read initial absorbance. Start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes. Determine the mean increase of absorbance per minute (dA/min).

### Calculation for sample start:

340 nm	dA/min. x 16030 = Activity (U/l)
Hg 334 nm	dA/min. x 16345 = Activity (U/l)

### Manual procedure for substrate start:

Wavelength	340nm (Hg 334nm)
Temperature	+37°C
Cuvette	1 cm light path
Zero adjustment	Air or distilled water

Serum/Plasma

R1	1000 $\mu$ l
Sample	10 $\mu$ l

Mix and incubate 1 to 5min. Then add:

R2	250 $\mu$ l
----	-------------

Mix, incubate for 60 sec., read initial absorbance. Start stopwatch simultaneously. Read again after exactly 1, 2 and 3 minutes. Determine the mean increase of absorbance per minute (dA/min).

### Calculation for substrate start:

340 nm	dA/min. x 20000 = Activity (U/l)
Hg 334 nm	dA/min. x 20390 = Activity (U/l)

## Measuring /reportable range:

Measuring range: up to 1200 U/l

Determine samples having higher activities via the rerun function. On instruments without rerun function, manually dilute such samples using 0.9% NaCl or distilled/deionised water (e.g. 1 + 9). Multiply the result by the appropriate dilution factor (e.g. 10)

## Expected values:

Adults: 72-182 U/l at 37°C

Quotient HBDH/LDH 0.63-0.81

A damage occurs at:

Myocardial lesion	HBDH/LDH	> 0.9
Kidney disease	HBDH/LDH	< 0.6

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes,  $\alpha$ -HBDH results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

## Analytical sensitivity (lower detection limit)

Detection limit: 3 U/l

The lower detection limit represents the lowest measurable  $\alpha$ -HBDH activity that can be distinguished from zero.

# Fluitest® $\alpha$ -HBDH

$\alpha$ -HYDROXYBUTYRATE DEHYDROGENASE

## **Imprecision:**

### **Serum**

Reproducibility was determined using controls and the following results were obtained (n=20):

Sample	With in run		
	Mean U/l	SD U/l	CV %
Sample 1	100	2.21	2.20
Sample 2	174	2.97	1.71
Sample 3	388	4.20	1.08

Sample	Run to run		
	Mean U/l	SD U/l	CV %
Sample 1	97.8	2.20	2.25
Sample 2	177	2.01	1.14
Sample 3	386	6.96	1.80

## **Method comparison:**

A comparison of the Analyticon Fluitest  $\alpha$ -HBDH (y) with a commercial obtainable assay (x) gave the following result (U/l):

$$y = 1.00x - 1.00; \quad r = 0.999$$

## **Quality control:**

Human Control Serum

Contronorm® Plus	5 x 5 ml	#1205
	20 x 5 ml	#1220
Controptath® Plus	5 x 5 ml	#1305
	20 x 5 ml	#1320

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

## **Calibration:**

S1: 0.9% NaCl

S2: Bio Cal® E                      10 x 3 ml                      #1430

## **Calibration frequency:**

Two-point calibration is recommended:

- after reagent lot change
- as required following quality control procedures

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