


CombiScan 100 Urine Analyzer



**Near patient
testing –
at it's best**




- CCD Colour Image Sensor
- LCD Display
- Throughput of 120 tests per hour
- Interface ports for connection to host computer, barcode reader and PC keyboard
- Urine colour compensation



The changeable strip holder can easily be removed for cleaning.



Simple changing of the thermopaper for the built-in printer. Rollwidth: 57 mm, diameter: 50 mm.



Connections on the back side of the Combi-Scan 100: ① = 7.5 V network connection, ② = USB connection, ③ = Keyboard

Features & Functions:

- Urine colour compensation
- Enter data via keyboard or barcode reader
- Variable format selection for printout of test results and time/data

Technical data:

- Type: CCD Colour Image Sensor
- Throughput: 120 strips/hour in fast mode, 40 strips/hour in convenient mode
- Display: 2 x 24 character LCD
Languages: English, German, French, Italian, Spanish
- Output options:
Languages: Chinese, English, German Spanish, French, Italian
Units: SI, Conventional, Arbitrary
- Printer: Low noise, built-in using 57 mm width and 50 mm diameter thermotransfer paper
- Memory: Last 500 results
- Interfaces: USB and PS2
- External Power Adapter: 7,5 V DC/3A

Optionals:

- Keyboard
- Barcode Reader

(All Modifications in view of technical progress are reserved)



Urine Strip Analyzer

Combi-Scan® 100

Art.-#

A93 009

Urine Test Strips

Combi-Screen® 11SYS

Combi-Screen® 11SYS PLUS

Pack size

100 tests

100 tests

Art.-#

93 100

94 100

Parameters:

Bilirubin, Urobilinogen, Ketones, Ascorbic Acid, Glucose, Protein, Blood, pH, Nitrite, Leukocytes, Specific Gravity

Combi-Screen® and Combi-Scan® are registered trademarks of Analyticon Biotechnologies AG.

Distributor information

**Analyticon
Biotechnologies AG**

Am Muehlenberg 10
35104 Lichtenfels - Germany

Fon +49 (0) 6454 7991-0
Fax +49 (0) 6454 7991-30

info@analyticon-diagnostics.com
www.analyticon-diagnostics.com