

# CombiScreen® Urine Control



## CombiScreen® Dip Check

Art.-No. 93010

Ready-to-use dipper control

2 x 15 ml (Level 1 + Level 2)

Open vial stability (at 2–8 °C) of 75 days or 20 determinations (whichever occurs first)

## CombiScreen® Drop Check

Art.-No. 93015

Ready-to-use dropper control

2 x 5 ml (Level 1 + Level 2)

Open vial stability of 18 months (2–8 °C) or 30 days (20–25 °C)

### Features:

- Based on human standard material
- Suitable for use in POC testing
- Shatter proof vials (polystyrene)
- 2D-barcode on tube vial for direct entry of lot number and target values into Urilyzer® 100 Pro and Urilyzer® 500 Pro instruments
- Target values for CombiScreen® test strips as well as other test strip brands
- Parameters: Bilirubin, Urobilinogen, Ketones, Glucose, Protein, Blood, pH, Nitrite, Leukocytes, Specific Gravity, Creatinine and hCG
- Qualitative hCG values

CombiScreen® urine controls cover a wide range of analytes, including pregnancy markers and Microalbumin. The controls are designed for use in manual and automated methods, to monitor the performance of a variety of urine test strips.

### Analyticon Biotechnologies AG

Am Muehlenberg 10  
35104 Lichtenfels - Germany  
Phone: +49 6454 7991-0  
info@analyticon-diagnostics.com  
www.analyticon-diagnostics.com

agile - affordable - accurate