

CombiScreen® Drop Check:

Compact and efficient solution
for your quality management

Ready-to-use urine control

Wide range of analytes including
pregnancy markers and microalbumin

2 x 5 ml (Level 1 + Level 2)

Open vial stability 18 month
(at 2-8°C) or 30 days (at 20-25°C)



Features:

- Based on human urine
- Parameters: Bilirubin, Urobilinogen, Ascorbic acid, Ketone, Glucose, Protein, Blood, pH, Nitrite, Leukocytes, Specific gravity, Microalbumin, Creatinine and hCG
- Target values for CombiScreen® urine test strips as well as other brands
- 2D barcode allows direct entry of batch number and target values during evaluation of the test strips on Urilyzer® 100 Pro and Urilyzer® 500 Pro devices
- Suitable for near-patient testing (POCT)



Analyticon Biotechnologies AG

www.analyticon-diagnostics.com | info@analyticon-diagnostics.com | +49 6454 7991 0