



CombiScreen Dip Check / CombiScreen Drop Check available now!

- Based on human standard material
- Suitable for use in POC testing
- Shatter proof vials (polystyrene)
- 2D-barcode on tube vial for direct entry of lot number and target values into Urilyzer[®] 100 Pro and Urilyzer[®] 500 Pro instruments
- Target values for CombiScreen[®] test strips as well as other test strip brands
- Parameters: Bilirubin, Urobilinogen, Ketones, Glucose, Protein, Blood, pH, Nitrite, Leukocytes, Specific Gravity, Creatinine and hCG
- Qualitative hCG values

[More Information here.](#)